

Clot Line

Coagulation Analyzers

CLOT SP

Semiautomatic coagulometer

Automatic clot detection

Programmed methods for PT, APTT, Fibrinogen

Thermostatized block for incubation of samples

and reagents

1 reading channel

Internal printer



CLOT 2B

Semiautomatic coagulometer

Automatic clot detection

Programmed methods for PT, APTT,

Fibrinogen

Thermostatized block for incubation of

samples and reagents

2 reading channels

Internal printer



CLOTTY

Semiautomatic coagulometer
Automatic clot detection
Programmed methods for PT, APTT, Fibrinogen
Thermostated block for incubation of
samples and reagents
1 reading channel



Reagents

Hemoplastin L - Quick Test PT
Hemos PTT - Ellagic Acid + Calcium Chloride
Hemofibrin - Bovine Thrombin + Imidazole Buffer
Hemotrol N
Hemotrol P
Coacal - Calibrator



RAL Técnica para el Laboratorio S.A. – 08970 Sant Joan Despí (Barcelona) SPAIN
ral@ral-sa.com – Tel +34 934808047 – Fax +34 933730092
www.ral-sa.com

